

Compatibility of Treff Pipette Tips with Pipettes

	Pipette tips (µl)					Filter tips (µl)							
	0.1-20	0.5-20	1-50	2-200	5-300	50-1000	0.1-1	0.5-10	1-20	2-20	5-100	5-200	50-1000
Treff®													
Transferpette® S													
0.1 - 1 µl	●						●						
0.1 - 2.5 µl	●	●	●				●	●	●				
0.5 - 10 µl	●	●	●				●	●	●				
2 - 20 µl				●	●					●	●		
10 - 100 µl				●	●					●	●		
20 - 200 µl				●	●					●	●	●	
100 - 1000 µl						●							●
Eppendorf®													
Reference®													
0.1 - 2.5 µl	●	●						●					
0.5 - 10 µl	●	●	●					●	●				
2 - 20 µl (grey)	●	●							●				
2 - 20 µl (yellow)				●						●		●	
10 - 100 µl				●						●	●	●	
50 - 200 µl				●		●					●	●	
100 - 1000 µl						●							●
Reference® 2													
0.1 - 2.5 µl	●	●	●				●	●	●				
0.5 - 10 µl	●	●	●				●	●	●				
2 - 20 µl (grey)	●	●	●					●	●				
2 - 20 µl (yellow)				●									
10 - 100 µl				●	●					●	●	●	
20 - 200 µl				●						●	●	●	
30 - 300 µl				●	●						●	●	
100 - 1000 µl						●							●
Research®													
0.1 - 2.5 µl	●	●						●					
0.5 - 10 µl	●	●						●					
2 - 20 µl (yellow)				●						●	●	●	
10 - 100 µl				●						●	●	●	
20 - 200 µl				●							●	●	
100 - 1000 µl						●							●
Research® plus													
0.1 - 2.5 µl	●	●						●					
0.5 - 10 µl	●	●	●					●	●				
2 - 20 µl (grey)	●	●	●						●				
2 - 20 µl (yellow)				●						●	●	●	
10 - 100 µl				●						●	●	●	
20 - 200 µl				●							●	●	
30 - 300 µl					●							●	
100 - 1000 µl						●							●
Xplorer®													
0.5 - 10 µl	●	●	●					●	●				
5 - 100 µl				●							●		
15 - 300 µl													
50 - 1000 µl						●							
Xplorer® plus													
0.5 - 10 µl	●	●	●					●	●				
5 - 100 µl				●							●		
15 - 300 µl													
50 - 1000 µl						●							●
BIOHIT®													
mLINE													
0.1 - 3 µl	●	●	●					●	●				
0.5 - 10 µl	●	●						●					
2 - 20 µl				●					●				
10 - 100 µl				●						●	●		
20 - 200 µl				●									
100 - 1000 µl						●							●

● = compatible
 □ = no information

The informations are based on general experiences and do not imply any warranty. The customer should always check the compatibility under his working conditons.

Status as of 2017

Compatibility of Treff Pipette Tips with Pipettes

	Pipette tips (µl)					Filter tips (µl)								
	0.1-20	0.5-20	1-50	2-200	5-300	50-1000	0.1-1	0.5-10	1-20	2-20	5-100	5-200	50-1000	
GILSON®														
Pipetman P														
0.2 - 2 µl	•	•	•					•	•					
1 - 10 µl	•	•	•					•	•					
2 - 20 µl				•						•				
20 - 100 µl				•							•	•		
50 - 200 µl				•								•		
200 - 1000 µl						•							•	
Pipetman Neo														
0.2 - 2 µl	•	•	•					•						
1 - 10 µl	•	•	•											
2 - 20 µl				•						•				
10 - 100 µl				•							•			
20 - 200 µl				•								•		
100 - 1000 µl						•							•	
Thermo														
FINNPIPETTE® Digital														
0.2 - 2 µl (micro)														
0.5 - 10 µl (micro)														
0.5 - 10 µl	•	•						•	•					
2 - 20 µl														
5 - 50 µl				•							•	•		
10 - 100 µl				•								•		
20 - 200 µl				•								•		
100 - 1000 µl						•							•	
FINNPIPETTE® Focus														
0.3 - 3 µl	•	•	•											
0.5 - 5 µl														
1 - 10 µl (pink)	•	•						•		•				
1 - 10 µl (yellow)				•						•				
FINNPIPETTE® Focus														
3 - 30 µl (turquoise)														
3 - 30 µl (yellow)														
5 - 50 µl (turquoise)														
5 - 50 µl (yellow)														
10 - 100 µl				•							•	•		
30 - 300 µl					•									
100 - 1000 µl						•							•	
Abimed														
Discovery														
0.5 - 10 µl	•	•						•	•					
2 - 20 µl				•						•				
10 - 100 µl				•										
20 - 200 µl				•								•		
100 - 1000 µl						•							•	
Labmate														
0.5 - 10 µl	•	•	•					•	•					
2 - 20 µl				•						•				
10 - 100 µl				•							•			
20 - 200 µl				•								•		
100 - 1000 µl						•							•	

= compatible
 = no information

The informations are based on general experiences and do not imply any warranty. The customer should always check the compatibility under his working conditons.

Status as of 2017